




## Clinical Trial Results




▶ Mean residual leukocyte count of 0.032 and 0.245 for BC and PRP method respectively



▶ Maximum hemolysis of 0.27% at end of storage



▶ Mean filtration time of 35 minutes



▶ Mean hemoglobin recovery of 49 g

## Specifications

USAGE CONDITIONS	
Whole blood collection volume	450ml ± 50ml (excluding CPD: 63ml)
Anticoagulant	CPD
Filtration temperature (during filtration)	Room temperature (~ 22° C)
Filtration method	Gravity
PERFORMANCE	
Residual WBCs	< 1 x 10 <sup>6</sup> / bag with conformance 90 %
RBC recovery	≥ 85 % (average)
Residual hemoglobin	≥ 40 g / bag
Hemolysis	< 0.8 % at end of storage
Filtration time	After holding at normal condition average 35 min

THE BRIDGE TO GOOD HEALTHCARE



Excellent and reliable

## Red Cell Filter System

BLOOD BAGS WITH IN-LINE LEUKOCYTE FILTER

For more information, please contact us:



Bionic Medizintechnik GmbH – a JMS corporation  
 Max-Planck-Str. 21 · D-61381 Friedrichsdorf, Germany  
 Telefon: +49(0)6172-75760 · Fax: +49(0)6172-757610  
 E-mail: info@bionic-jms.com · www.bionic-jms.com

Rights for technical modifications reserved  
 JMS-RCFE01 2017-04

# Red Cell In-Line Filtration System

## Advantages at a glance:

- ▶ Global In-Line leukocyte filter for red cell components with excellent performance
- ▶ Standard procedures for collection, filtration and storage
- ▶ First class filter media technology proven in markets

## In-Line Filter

Premium In-Line leukocyte filter for red cell components with soft housing


## Filtration Procedure

Standard filtration procedure and method for easier handling

Sepacell is a trademark of Asahi Kasei Medical Co., Ltd.




# Blood Components Separation Procedure



■ Red Blood Cells (RBC)    
 ■ Plasma    
 ■ Buffy Coat    
 ■ Platelet Concentrate (PC)

# Filtration Procedure



**Filtration:** Turn Red Cells bag and place at height of 1.4 m. Mix well and open clamps to commence filtration

**Filtration:** Close red clamp to cease filtration once flow has stopped

**Completion:** Discard filtration portion