


Don't just seal, SURESEAL



CENTRON

SureSeal™ SE260

TUBE SEALER

The professional choice for heavy-duty tube sealing

SureSeal™ SE260 is AC powered benchtop tube sealer with integrated sealing head and optional sealing hand unit. New sealing technology and enforced reliability make SE260 ideal for most demanding environments.

■ Features


- Premium sealing quality
- Simple operation, fast sealing, no warm-up required
- Sealing hand unit ready
- Automatic tubing detection
- LED indicator to show sealing progress
- Applicable to wide range of blood bag tubing
- Head protector / Splash guard for user safety

Centron

Smart Solutions for Smarter Bloodbanking

New state-of-the-art sealing head

- New sealing technology ensures safe sealing of even thick and sensitive tubings with wide sealed surface.
- With automatic tubing detection, sealing will automatically start when tubing is placed in the sealing head.
- Semitransparent head protector for safe and easy operation (detachable, easy to clean)
- Applicable tubing diameter 2-6mm (1/4")




Splashguard

- protects users from accidental blood splash
- Detachable, easy to clean

Indication lamp

- shows sealing progress




Wide sealed surface with notch at the center

- Reliable sealing quality prevents leakage
- Easy to segment



Sealing hand unit ready

- Additional plug-in type *SureSeal™* sealing hand unit
- Ergonomic design with standard 1.8m (6ft) cable
- Satisfying various operating conditions



Specification


Product	TUBE SEALER
Model	SE260
Power	100-240V~ 50/60Hz 250W
Dimension	123(W)mm/302(L)mm/243(H)mm

Information



PENDING

DISTRIBUTOR



Centron Centron Corporation

Address: 44 Beolmal-ro, Dongan-gu, Anyang, Gyeonggi-do,
Postal Code 431-804 / Republic of Korea
Phone: +82-31-422-7890 / **Fax:** +82-31-422-7866
E-mail: info@centrons.com
Website: www.centrons.com



※ The specification and design are subject to change without prior notice for product upgrades.