



Mobile heat for
infusions and textiles

The Barkey warming center II Adjustable heat for infusions and textiles according to your requirements

- Adjustable temperature range
- High energy efficiency
- User-friendly mobility
- Battery-operated display (optional)
- Large capacity
- Optional accessories
- Simple to use
- Overtemperature alarm

warming center II



Barkey GmbH & Co. KG

Gewerbestrasse 8
33818 Leopoldshoehe
Germany
www.barkey.de

Temperature range directly adjustable by the user



from +35° to +40°C for infusion fluids



and from +41° to +55°C for textiles

High energy efficiency through intelligent heating technology and energy-saving insulation – leads to a 50% time saving when heating fluids in comparison to the previous model. 40 infusion bottles (1000 ml bottles at a starting temperature of 21°C) are heated to 37°C within approx. 8 hours.

User-friendly mobility thanks to a large transport handle on the rear and four stable castors (2 with locking brakes).

Universally applicable
For heating infusions and textiles.



Use with textiles

Large capacity – each of the 2 stable drawers provides space for up to 20 oval 1000 ml or 18 round 1000 ml bottles and can take a load of up to 22 kg.

Mechanical blocking of drawer pull-out for safe handling.

Overtemperature alarm both static and dynamic is also part of the standard equipment as well as an alarm for when a drawer is left open for more than 1 minute.

Optional accessories are also available. Increase the functionality of your warming center II with a data logger and/or a partitioning system to safeguard bottles against tipping over.

Battery-operation of the display possible with a battery module available separately.

Type	Barkey warming center II
Power supply	230 V 50/60 Hz or 100-120 V 50/60 Hz
Standby	< 1,5 W
Temperature setting	selectable from +35° to +55°C
Default setpoint	+39° C
Battery type (optional)	NIMH, 6 VDC, 2.100 mAh,
Typical capacity	40 bottles each with 1000 ml,
Vol.	approx. 80 litres
Weight empty	59 kg
External dimensions	529 x 912 x 698 mm (W x H x D)
Drawer size (internal)	365 x 275 x 410 mm (W x H x D)
Drawers	Aluminium/stainless steel with perforated base, soft-close
Load capacity per drawer	max. 22 kg
Protection Class	I
MDD classification	II a
Certification	CE 0123

Intended purpose

The warming center II is used for warming up fluids intended for perfusion contained in bottles (up to 1000 ml) or bags (up to 5000 ml), textiles (bedsheets, etc.) and small items (syringes, instruments, etc.) to a temperature of + 35°C to + 55°C and maintaining their heat.

Sold by:

www.barkey.de